

PMS Parts Reference Guide

Part Number: SP00218

Rev: -

Date: 31 May 2001

© 2001 Draeger Medical, Inc.

Incubator 8000


PMS Parts Reference Guide

This document is intended as a reference guide to assist in the replacement of routine maintenance parts.


Incubator 8000

| <u>Description</u> | <u>Part Number</u> | <u>Comments</u> |
|--------------------|---|--|
| 6 month | <i>No Kit</i> 2M07839 | Oil 10 drops |
| 24 month | 4114770 8402926 2M19622 1828142 8301856 8409401 | Filter 1x Filter 1x D-Ram 1x Storage Battery 1x DISS Adapter (IC only) 1x |
| 72 month | <i>No kit</i> 8400496 8402745 <i>Or</i> 8406678 <i>Or</i> 8411142 | Spring (front flap catch) 2x Pressure Regulator 1x Reg Repair Set (Knocks) Reg Repair Set (Lorch) |


Incubator




Motor




Plug – In Unit




Power Pack PCB




CPU PCB



Regulator



Front Flap



DrägerService®

Dräger

DrägerService is a division of
Draeger Medical, Inc.
3122 Commerce Drive
Telford, PA 18969
Tel: (215) 721-5402
(800) 543-5047
Fax: (215) 721-5784
Web: www.draegermedical.com
Printed in the U.S.A.